



## MRI SAFETY INFORMATION

### OPTIMIZER Smart Mini Implantable Pulse Generator

	MR - Conditional
---	------------------

**MR scans are permissible at 1.5 T and 3 T with local RF transmit-receive head and extremity coils only. The use of the RF Body coil can result in serious injury to the patient.**

The OPTIMIZER Smart Mini IPG is Magnetic Resonance (MR) Conditional, and patients with this device may be scanned safely with magnetic resonance imaging (MRI) **if all the requirements for the implanted components and for scanning are met.**

### Restrictions for the Patient and the Implanted System

- The OPTIMIZER Smart Mini IPG must be implanted with two ventricular leads (and one optional atrial lead) of the following models:

Model	Manufacturer	Lengths (cm)
CapSureFix Novus MRI™ SureScan™ 5076	Medtronic	52, 58
CapSureFix Novus MRI™ SureScan™ 5086	Medtronic	52, 58
SelectSecure MRI SureScan 3830	Medtronic	59
2088TC Tendril STS	Abbott (St Jude Medical)	52, 58
Solia S	Biotronik	53, 60
Ingevity 7741 and 7742	Boston Scientific	52, 59
Ingevity + 7841 and 7842	Boston Scientific	52, 59

When combined, the OPTIMIZER Smart Mini IPG and such leads constitute an MR conditional device system.

**WARNING:** Only the lead lengths above have been evaluated for MRI compatibility and individual scan parameters. The system is considered MR-Conditional only if used with the listed leads and lengths.

- There are no other active or abandoned cardiac implants (e.g., lead extensions, lead adapters or abandoned leads) in the patient's body.

**WARNING:** Do not bring any system components that are not marked MR-safe or MR-conditional into the MRI suite.

- At least six (6) weeks have elapsed since the OPTIMIZER Smart Mini IPG and/or lead implantation and/or any electrode revision or surgical modification.

- The device system is implanted pectorally.
- The OPTIMIZER Smart Mini IPG is programmed to OOO mode before the MRI scan.
- Patient does not have elevated body temperature or compromised thermoregulation at time of scan.

**WARNING:** Do not scan a patient with an elevated body temperature.

## Requirements of the MRI Scanner

- Use of a clinical, hydrogen-atom MRI scanner with horizontal cylindrical closed-bore magnet, and a static magnetic field strength of **1.5 or 3 Tesla**.
- **Use only transmit-receive head and extremities local RF coils** that are not positioned over the implant location.
- Maximum spatial gradient of the magnetic field of 40 T/m or 4000 Gauss/cm.
- The slew rate of the MRI scanner's gradient fields must not exceed 200 T/m/s per axis.

**WARNING:** scanning under other conditions may result in severe patient injury, death, or device malfunction.

## Restrictions During the MRI scan

- The specific absorption rate must not exceed limits for local transmit-receive coils as defined in IEC 60601-2-33, for example 3.2 W/kg for the head.
- Emergency equipment for resuscitation must be kept at hand and properly certified staff must be available.
- Patient must be continuously monitored by pulse oximetry and electrocardiography (ECG).

## Image Artifacts

In non-clinical testing, the maximum image artifact size was seen on the gradient echo pulse sequence at 1.5 T and extends by approximately 6.7 cm from the boundary of the implant.

In non-clinical testing, the maximum image artifact size was seen on the gradient echo pulse sequence at 3 T and extends by approximately 4.6 cm from the boundary of the implant.

**Impulse Dynamics (USA), Inc.**  
50 Lake Center Executive Parkway  
401 Route 73 N, Building 50, Suite 100  
Marlton, NJ 08053-3449  
(856) 642-9933 [www.impulse-dynamics.com](http://www.impulse-dynamics.com)

