



**Vesta Charger for
OPTIMIZER™ Smart Mini
System**

INSTRUCTIONS FOR USE



Impulse Dynamics (USA) Inc.
Suite 100
50 Lake Center Executive Parkway
401 Route 73 N Bldg. 50
Marlton, NJ 08053-3425



Impulse Dynamics Germany GmbH
MAC Main Airport Center
Unterschweinstiege 2-14
60549 Frankfurt am Main
Germany

OPTIMIZER™ and CCM™ are trademarks of Impulse Dynamics.

OPTIMIZER® is a US registered trademark property of Impulse Dynamics.

The Vesta Charger conforms to the essential requirements of the Radio Equipment Directive (RED) 2014/53/EU.

The information provided in this document may change without prior notice.

No part of this manual may be reproduced or transmitted in any form or by any method, including electronic and mechanical means, without prior written express consent from Impulse Dynamics.

The OPTIMIZER Smart Mini system and the CCM technology are protected by several U.S. Patents. For an up-to-date list of relevant patents and patent applications, visit our patents page:

<http://www.impulse-dynamics.com/us/patents>

Please read the documentation provided completely before you use the device.

NOTICE: Any serious incident should be reported to the manufacturer Impulse Dynamics by emailing QualityComplaints@impulse-dynamics.com. According to MDR 2017/745, a 'serious incident' means any incident that directly or indirectly led, might have led, or might lead to any of the following:

- a) The death of a patient, user, or other person,
- b) The temporary or permanent serious deterioration of a patient's, users, or other person's state of health. A serious deterioration in the health of the subject, that resulted in any of the following:
 - i. Life-threatening illness or injury,
 - ii. Permanent impairment of a body structure or a body function,
 - iii. Hospitalization or prolongation of patient hospitalization,
 - iv. Medical or surgical intervention to prevent life-threatening illness or injury or permanent impairment to a body structure or body function,
 - v. Chronic disease,
- c) A serious public health threat. A public health threat is an event which could result in imminent risk of death, serious deterioration in a person's state of health, or serious illness, that may require prompt remedial action, and that may cause significant morbidity or mortality in humans, or that is unusual or unexpected for the given place and time.


















Revision 00, Date of Issue: 2022-07-15

TABLE OF CONTENTS

EXPLANATION OF SYMBOLS ON LABELS	I
1.0 THE VESTA CHARGER SYSTEM	1
1.1 Description	1
1.2 User Profile and Training	1
1.3 Charging Method	2
1.4 System Components	2
1.5 Features	3
1.6 Overview of the Screens Displayed by the Vesta Charger	3
1.6.1 Screens Displayed When Connected to the AC Adapter	3
1.6.2 Screens Displayed When Pairing with the OPTIMIZER Smart Mini IPG	5
1.6.3 Screens Displayed When Charging the OPTIMIZER Smart Mini IPG	5
1.6.4 Screens Displayed After the Detection of an Alert Condition	8
1.6.5 Info Screens	10
1.7 Pairing the Vesta Charger with the OPTIMIZER Smart Mini IPG	11
1.8 Charging the Vesta Charger	12
1.9 Charging the OPTIMIZER Smart Mini IPG	13
1.9.1 Early Termination of Charging Session	15
1.10 Charging the OPTIMIZER Smart Mini IPG in Special Charge Mode	15
1.11 Vesta Charger Placement When Not Being Used for Device Charging	16
1.12 Frequency of Charging Sessions	16
1.13 Communications	17
1.13.1 Communications with the OPTIMIZER Smart Mini IPG	17
1.14 Call Doctor Alert Codes	17
1.14.1 Call Doctor Alert Code Attributes	18
1.14.2 Call Doctor Alert Code Definitions	18
1.15 FCE Charger	19
1.15.1 Screens Displayed When FCE Charger is Connected to AC Adapter	19
1.15.2 Charging the OPTIMIZER Smart Mini IPG using the FCE Charger	19
1.15.3 Charging the OPTIMIZER Smart Mini IPG in Unpaired Charge Mode	21
1.16 Cleaning	22
1.17 Maintenance	23
1.18 Storage and Handling	23
1.19 Disposal	23
APPENDIX I	25
Electromagnetic Immunity	25
Electromagnetic Immunity of the Vesta Charger	25
Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the Vesta Charger	27
Electromagnetic Emissions	28

Electromagnetic Emissions from the Vesta Charger	28
APPENDIX II	31
Wireless Technology	31
Vesta Charger Wireless Nominal Specifications	31
Quality of Service (QoS) for Wireless Technology	32
Wireless Security Measures	33
Troubleshooting for Wireless Coexistence Issues	33

EXPLANATION OF SYMBOLS ON LABELS

Symbol	Description
	CE Conformity marking, 0344 - Notified Body Number
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult instructions for use
	Do not use if package is damaged
	Storage and transport temperature limits
	Date of manufacture
	Manufacturer
	Authorized representative in the European Community
	Catalogue number
	Serial number
	Refer to instruction manual/booklet
	Caution, consult instructions for use
	Item not to be disposed via the municipal waste collection system of any member state of the European Union
	Type BF applied part
	Non-ionizing electromagnetic radiation
IP22	Protected against the ingress of solid foreign objects greater than 12.5 mm (0.5 in) in width Protected against the ingress of vertically falling water drops when the enclosure is tilted at an angle of 15° from its normal position

THIS PAGE INTENTIONALLY LEFT BLANK

1.0 THE VESTA CHARGER SYSTEM

1.1 Description

The Vesta Charger is designed to charge the battery of the OPTIMIZER Smart Mini IPG with only minimal patient intervention while ensuring patient safety and maintaining proper operation of the IPG during the charging process.

In addition, the Vesta Charger is programmed to display alerts and other messages that may require action by the patient (e.g., Call Doctor Alert Codes that require the patient to contact the doctor, reminders to charge your implanted device, etc.).

The Vesta Charger has a permanently attached charging wand and is powered by a rechargeable battery. To recharge this battery, the Vesta Charger System includes a Cell-Con AC Adapter (Input: 100–240 VAC, 50-60 Hz, 0.2 A; Output: 4.2 V, 1.3 A).

The Vesta Charger is a Class I, Type BF device, classified as ordinary equipment suitable for continuous operation, with short-time loading, within the patient environment.

Caution: The Vesta Charger is subject to interference from other electrical devices operated in the vicinity. Portable and mobile Radio Frequency (RF) equipment are especially prone to impair the normal function of the charger. If Vesta Charger is not operating as expected, such interference has to be taken into account.

The Vesta Charger communicates with the OPTIMIZER Smart Mini IPG at a frequency range of 402 MHz to 405 MHz (MedRadio frequency band). The communication range of the Vesta Charger is between zero and at least 1.5 m (5 ft).

The Vesta Charger charges the OPTIMIZER Smart Mini IPG at a frequency range of 13.56 MHz.

When the distance between the Charging Wand and the OPTIMIZER Smart Mini IPG is between 0.5 cm and 2.0 cm, the Vesta Charger should be able to recharge the OPTIMIZER Smart Mini IPG with a battery charge of 10% to 90% battery charge in less than 2.5 hours with the Vesta Charger's charging current fixed at 90 mA \pm 10%.

When the distance between the Charging Wand and the OPTIMIZER Smart Mini IPG is > 2.0 cm, the Vesta Charger should be able to recharge the OPTIMIZER Smart Mini IPG in less than 4 hours with starting and ending IPG battery charge levels shown in the **Table 1**. In such an instance, the Vesta Charger's battery may become depleted before the rechargeable battery in the OPTIMIZER Smart Mini IPG is fully charged

Table 1: IPG Charge Levels Attained with Fully Charged Vesta Charger

Charging Wand and IPG Distance	Starting IPG Battery Charge Level	Ending IPG Battery Charge Level
> 2.0 cm, \leq 3.5 cm	10%	80%
> 3.5 cm, \leq 4.0 cm	10%	70%

When fully charged, the Vesta Charger should be able to perform two IPG charging cycles, charging the IPG battery from 10% to 90% each time, when the distance between the Charging Wand and the OPTIMIZER Smart Mini IPG is between 0.5 cm and 2.0 cm.

1.2 User Profile and Training

The operators of the Vesta Charger system include patients, physicians (and the trained medical personnel who assist them), and Impulse Dynamics representatives. Physicians, medical personnel, and Company representatives who operate the Vesta Charger system should be familiar with the operation of electronic medical equipment, particularly the operation of implanted medical devices.

Physicians and medical personnel can participate in a Company-sponsored training program that will provide theoretical and hands-on training regarding the technology, device features, and detailed operating instructions for the Vesta Charger.

Patient training on the use of the Vesta Charger will be provided by Impulse Dynamics Representatives post-implant.

1.3 Charging Method

The charging method utilized by the Vesta Charger to charge the battery of the OPTIMIZER Smart Mini IPG is called inductive energy transfer. Since magnetic fields can penetrate human tissues with nearly no attenuation, inductive energy transfer is the only practical transcutaneous recharging method.

The manner in which inductive energy transfer is used to charge the battery of the OPTIMIZER Smart Mini IPG is as follows:

1. Electrical energy from the battery of the Vesta Charger passes through a primary coil connected to the electronic circuitry of the charger that converts it into an oscillating electromagnetic field.
2. When a primary coil is placed in close proximity to a secondary coil, the oscillating electromagnetic field generated by a primary coil is picked up by a secondary coil.
3. The secondary coil that picks up the oscillating electromagnetic field is connected to the electronic circuitry of the implant that converts it back into electrical energy. That electrical energy is used to charge the battery of the OPTIMIZER Smart Mini IPG.

1.4 System Components

The Vesta Charger System consists of the following components:



Figure 1: Vesta Charger System Components

- **Vesta Charger** (with attached charging wand and charging wand cable clip) – used to charge the OPTIMIZER Smart Mini IPG.
- **AC Adapter** – used to charge the internal battery of the Vesta Charger.
- **EU/US Plug Adapters** – plug adapters for the AC Adapter, allowing the AC Adapter to be connected to wall outlets in the EU and US.
- **Carrying Case** – used to store and transport the Vesta Charger System.

1.5 Features

The Vesta Charger has the following features:

- **Graphical Display:** Display screen used by the Vesta Charger to communicate information to the patient
- **Power Button:** Press-button switch used to initiate and terminate charging of the OPTIMIZER Smart Mini IPG and to silence alerts displayed by the Vesta Charger
- **Buzzer:** An internal buzzer that produces beeping tones to inform the patient of a condition that requires action
- **Charging Wand:** Wand containing a coil and circuitry used by the Vesta Charger for charging as well as short-range communications with the OPTIMIZER Smart Mini IPG
- **Radio Transceiver:** Device used by the Vesta Charger for long-range communications [between zero and at least 1.5 m (5 ft)] with the OPTIMIZER Smart Mini IPG

1.6 Overview of the Screens Displayed by the Vesta Charger

The Vesta Charger displays a different screen for each operational state. This section presents an overview of each screen displayed by the Vesta Charger.

1.6.1 Screens Displayed When Connected to the AC Adapter

1.6.1.1 Charger Self-Charge Status Screen

This screen is displayed whenever the AC Adapter is connected to the Vesta Charger. The number of bars shown on the battery icon will vary depending on the current level of charge in the Vesta Charger battery (see **Table 2**).

Table 2: Vesta Charger Battery Charge Levels

Charger Battery Icon (When Not Charging or Charge Complete)	Charger Battery Icon (When Charging)	Charger Battery Charge Level
1 bar	1 flashing bar	Below 25%
2 bars	2 bars, last one flashing	Between 25% and 50%
3 bars	3 bars, last one flashing	Between 50% and 75%
4 bars	4 bars, last one flashing	Above 75%



Figure 2: Charger Self-Charge Status Screen

1.6.1.2 Charger Self-Charge Success Screen

This screen is displayed either when the AC Adapter has successfully completed charging the internal battery of the Vesta Charger, when the AC Adapter is connected to the Vesta Charger and the battery charge level of the Vesta Charger is above 75%, or when the AC Adapter is charging the Vesta Charger and the AC Adapter current is less than 50 mA.




Figure 3: Charger Self-Charge Success Screen

1.6.1.3 IPG Data Download Screen

This screen is displayed whenever the Vesta Charger is actively attempting to download data from the OPTIMIZER Smart Mini IPG. The encrypted data downloaded from the device includes information regarding the current status of your IPG, statistical information regarding its operation, and any active alerts that require action.

This is the first screen displayed after the AC Adapter is connected to the Vesta Charger and then plugged into the wall outlet.




Figure 4: IPG Data Download Screen

1.6.1.4 IPG Data Download Success Screen

This screen is displayed whenever the Vesta Charger has successfully completed downloading data from the OPTIMIZER Smart Mini IPG.

This is the second screen displayed after the AC Adapter is connected to the Vesta Charger and then plugged into the wall outlet.



Figure 5: IPG Data Download Success Screen

1.6.1.5 IPG Data Download Error Screen

This screen is displayed whenever the Vesta Charger has not successfully completed downloading data from the OPTIMIZER Smart Mini IPG.



Figure 6: IPG Data Download Error Screen

1.6.2 Screens Displayed When Pairing with the OPTIMIZER Smart Mini IPG

1.6.2.1 Charger/IPG Pairing Screen

This screen is displayed whenever the Vesta Charger is actively pairing with the OPTIMIZER Smart Mini IPG.




Figure 7: Charger/IPG Pairing Screen

1.6.2.2 Charger/IPG Pairing Success Screen

This screen is displayed whenever the Vesta Charger has successfully paired with the OPTIMIZER Smart Mini IPG. The display of this screen is accompanied by 3 short beeping tones.



Figure 8: Charger/IPG Pairing Success Screen

1.6.2.3 Charger/IPG Pairing Error Screen

This screen is displayed whenever an error has occurred during the pairing of the Vesta Charger and the OPTIMIZER Smart Mini IPG.



Figure 9: Charger/IPG Pairing Error Screen

1.6.3 Screens Displayed When Charging the OPTIMIZER Smart Mini IPG

1.6.3.1 IPG Data Download Screen

This screen is displayed whenever the Vesta Charger is actively downloading data from the OPTIMIZER Smart Mini IPG.

This is the first screen displayed after pressing the Power Button on the Vesta Charger to begin a charging session.




Figure 10: IPG Data Download Screen

1.6.3.2 IPG Data Download Success Screen

This screen is displayed whenever the Vesta Charger has successfully completed downloading data from the OPTIMIZER Smart Mini IPG. The display of this screen is accompanied by 3 short beeping tones.

If the Vesta Charger has successfully completed downloading data from the OPTIMIZER Smart Mini IPG, this is the second screen that is displayed after pressing the Power Button on the Vesta Charger to begin a charging session.




Figure 11: IPG Data Download Success Screen

1.6.3.3 IPG Data Download Error Screen

This screen is displayed whenever the Vesta Charger has not successfully completed downloading data from the OPTIMIZER Smart Mini IPG. The display of this screen is accompanied by 3 long beeping tones.

If the Vesta Charger is unable to establish coupling with the OPTIMIZER Smart Mini IPG, this is the second screen that is displayed after pressing the Power Button on the Vesta Charger to begin a charging session.



Figure 12: IPG Data Download Error Screen

1.6.3.4 Charging IPG Status Screen

This screen is displayed whenever the Vesta Charger has successfully coupled with the OPTIMIZER Smart Mini IPG and is charging the implanted device.

If the Vesta Charger has successfully coupled with the OPTIMIZER Smart Mini IPG, this is the third screen displayed after pressing the Power Button on the Vesta Charger.

The number of bars shown on the Vesta Charger Battery icon (on the left) and the IPG Battery icon (on the right) will vary depending on the current level of charge in each battery (see **Tables 3 and 4**).

Table 3: Vesta Charger Battery Charge Levels

Vesta Charger Battery Icon	Charger Battery Charge Level
1 bar	Below 25%
2 bars	Between 25% and 50%
3 bars	Between 50% and 75%
4 bars	Above 75%

Table 4: OPTIMIZER Smart Mini IPG Battery Charge Levels

IPG Battery Icon	IPG Battery Charge Level
1 flashing bar	Below 25%
2 bars, last one flashing	Between 25% and 50%
3 bars, last one flashing	Between 50% and 75%
4 bars, last one flashing	Above 75%




Figure 13: Charging IPG Status Screen

1.6.3.5 Charging IPG Coupling Error Screen

This screen is displayed whenever the Vesta Charger is unable to establish coupling with the OPTIMIZER Smart Mini IPG. The display of this screen is accompanied by 3 long beeping tones.

If your Vesta Charger is unable to establish coupling with the OPTIMIZER Smart Mini IPG, this is the third screen that is displayed after pressing the Power Button on the Vesta Charger.




Figure 14: Charging IPG Coupling Error Screen

1.6.3.6 IPG Charging Successfully Completed Screen

This screen is displayed whenever the Vesta Charger has successfully completed charging the battery of the OPTIMIZER Smart Mini IPG.




Figure 15: IPG Charging Successfully Completed Screen

1.6.3.7 Charging IPG Timeout Error Screen

This screen is displayed by the Vesta Charger whenever the charging duration of the OPTIMIZER Smart Mini IPG exceeds 5 hours ± 5 minutes.



Figure 16: Charging IPG Timeout Error Screen

1.6.3.8 Charging IPG Temperature Error Screen

This screen is displayed by the Vesta Charger whenever one of the following conditions occurs:

- The reported temperature of the OPTIMIZER Smart Mini IPG at the beginning of the charging session is outside the accepted range.
- The charging session is suspended due to the temperature of the OPTIMIZER Smart Mini IPG remaining consistently high for more than 10 minutes.




Figure 17: Charging IPG Temperature Error Screen

1.6.3.9 Power Supply Error Screen

This screen is displayed whenever the AC Adapter is connected to the Vesta Charger while it is charging the OPTIMIZER Smart Mini IPG.




Figure 18: Power Supply Error Screen

1.6.3.10 Charge Session Cancellation Screen

This screen is displayed whenever the button on the Vesta Charger is pressed while it is charging the OPTIMIZER Smart Mini IPG. The display of this screen is accompanied by 3 short beeping tones.

This screen is displayed just before the Vesta Charger shuts off.




Figure 19: Charge Session Cancellation Screen

1.6.4 Screens Displayed After the Detection of an Alert Condition

1.6.4.1 Low Charger Battery Alert Screen

This screen is displayed whenever the Vesta Charger's battery charge level drops below 10%. The display of this screen is accompanied by short beeping tones.



Figure 20: Low Charger Battery Alert Screen

1.6.4.2 Long Time Without Charging IPG Alert Screen

This screen is displayed whenever the Patient Alert “Battery Recharge Reminder” is enabled using the Optimizer Smart Mini Programmer application and the number of days since the OPTIMIZER Smart Mini IPG was last charged has exceeded the number of days set for this Patient Alert. The display of this screen is accompanied by short beeping tones.



Figure 21: Long Time Without Charging IPG Alert Screen

1.6.4.3 Long Time Without Downloading Data From IPG Alert Screen

This screen is displayed whenever the Patient Alert “Long Time Without Communicating with the IPG” is enabled using the Optimizer Smart Mini Programmer application and the number of days since the last successful communication between the Vesta Charger and the OPTIMIZER Smart Mini IPG has exceeded the number of days set for this Patient Alert. The display of this screen is accompanied by short beeping tones.



Figure 22: Long Time Without Downloading Data From IPG Alert Screen

1.6.4.4 Abnormal Condition Error Screen

This screen is displayed whenever an abnormal condition is detected in the OPTIMIZER Smart Mini IPG or the Vesta Charger. The display of this screen is accompanied by 3 long beeping tones.



Figure 23: Abnormal Condition Error Screen

1.6.4.5 Call Doctor Alert Screen

This screen is displayed whenever a Call Doctor Patient Alert that is enabled by the Optimizer Smart Mini Programmer application has been activated. The letter displayed is specific to the model of the implanted IPG. The display of this screen is accompanied by short beeping tones.



Figure 24: Call Doctor Alert Screen

1.6.4.6 Snooze Buzzer Alert Screen

This screen instructs the patient to press the button on the Vesta Charger to silence the beeping tone associated with the activated alert.

It is the screen that is displayed after the alert screen of a newly activated alert.




Figure 25: Snooze Buzzer Alert Screen

1.6.4.7 Snooze Alert Screen

This screen instructs the patient to press the button on the Vesta Charger to snooze an alert.

This screen is displayed after the alert screen if the Vesta Charger is used outside the scheduled Patient Alert Delivery period set by the Optimizer Smart Mini Programmer application (usually between 08:00 and 21:00) or when an alert that was previously activated is retriggered.



Figure 26: Snooze Alert Screen

1.6.5 Info Screens

The Vesta Charger displays the Info Screens when the following conditions are met:

- The AC Adapter is connected to the Vesta Charger.
- The **Power Button** is continuously pressed until a beeping tone is heard and then released (usually more than 5 seconds and less than 10 seconds).

1.6.5.1 First Info Screen

When the **Power Button** is released, the First Info screen displays the following information:

- The list of active and snoozed Call Doctor Alert Codes
- The IPG model code
- The battery charge level of the IPG after the completion of its last charge session
- The date and time of the last successful charge of the IPG

Note: The date format is (DD/MM/YY) and the time format is 24 hours.




Figure 27: First Info Screen

1.6.5.2 Second Info Screen

After the display of the First Info screen, the Second Info screen displays the following information:

- The signal strength level during the last successful IPG data download session
- The date and time of the last successful IPG data download session

Note: The date format is (DD/MM/YY) and the time format is 24 hours.



Figure 28: Second Info Screen

1.7 Pairing the Vesta Charger with the OPTIMIZER Smart Mini IPG

The pairing of the Vesta Charger with the OPTIMIZER Smart Mini IPG ensures that the communication and charging information received by the Vesta Charger is securely encrypted and unique to a specific implanted device.

During the pairing process, the Vesta Charger uses short-range communication to search for a device to pair with and creates an encryption key once a compatible device model has been found. This encryption key is stored and used by the Vesta Charger for all its subsequent communications sessions with the paired device.

To pair the Vesta Charger with the OPTIMIZER Smart Mini IPG, perform the following steps:

1. Determine the location of the OPTIMIZER Smart Mini IPG (typically right upper chest area) and then place the charging wand directly over the OPTIMIZER Smart Mini implant site (over the patient's clothes).
2. Place a pairing magnet (or a standard pacemaker magnet) to the left of the **Power Button** on the Vesta Charger. **See Figure 29.**

Note: A Vesta Charger being used for the first time does not require the use of a magnet during the pairing process.




Figure 29: Pairing Magnet on Vesta Charger

3. Begin the pairing process by pressing the **Power Button**, holding the button down for 1-2 seconds, and then releasing it.

- The Charger/IPG Pairing screen is displayed while the Vesta Charger is actively attempting to pair with the OPTIMIZER Smart Mini IPG. **See Figure 30.**




Figure 30: Charger/IPG Pairing Screen

- When the pairing process has been completed, the Vesta Charger will emit 3 short beeping tones and display the Charger/IPG Pairing Success screen. **See Figure 31.**



Figure 31: Charger/IPG Pairing Success Screen

- Remove the pairing magnet from the Vesta Charger.

1.8 Charging the Vesta Charger

Note: When the Vesta Charger is not being used to charge their implanted device, advise patients to always keep it connected to its AC Adapter and the AC Adapter plugged into the wall outlet. This keeps the battery of the Vesta Charger fully charged and ready to be used the next time they need to charge their implanted OPTIMIZER Smart Mini IPG.

Note: Charging the Vesta Charger and charging the OPTIMIZER Smart Mini IPG CANNOT be done simultaneously. Always charge the internal battery of the Vesta Charger before attempting to charge the battery of the OPTIMIZER Smart Mini IPG.

Note: Inspect the AC Adapter for any damage before each use. Contact your Impulse Dynamics representative if a replacement AC Adapter is needed.

Warning: Only use the AC Adapter provided with the Vesta Charger to charge the battery in the Vesta Charger. Otherwise, damage to the Vesta Charger may result.

To connect the AC Adapter to the Vesta Charger and begin charging its internal battery, perform the following steps:

- Turn the Vesta Charger around so that the back of the charger is facing up.
- Remove the protective cover flap from the power input connector located next to the base of the charging wand cable.
- Obtain the AC Adapter from the Carrying Case and rotate its DC output connector until the red dot on its connector is visible.
- Line up the red dot on the DC output connector of the AC Adapter with the red line on the power input connector of the Vesta Charger and then insert in the DC output connector into the power input connector. **See Figure 32.**




Figure 32: Connection of DC Connectors

- Attach the location-specific Plug Adapter to the AC Adapter and then plug the AC Adapter into the wall outlet to begin charging the internal battery of the Vesta Charger.

When the Charging Self-Charge Success screen is displayed on the Vesta Charger screen, the battery in the Vesta Charger is fully charged. **See Figure 33.**




Figure 33: Charger Self-Charge Success Screen

To disconnect the AC Adapter from the Vesta Charger, perform the following steps:

1. Unplug the AC Adapter from the wall outlet.
2. Hold and pull back on the metal sleeve of the DC output connector to disconnect it from the Vesta Charger.
3. Replace the protective cover flap over the power input connector of the Vesta Charger.

1.9 Charging the OPTIMIZER Smart Mini IPG

Warning: Failure to recharge the OPTIMIZER Smart Mini IPG as required may cause it to shut down when the battery is depleted, suspending CCM therapy delivery.

Note: The Vesta Charger cannot be used to charge the OPTIMIZER Smart Mini IPG until the AC Adapter is disconnected from the Vesta Charger.

Note: The Vesta Charger should not be operated close to other electronic equipment. If sufficient spatial separation cannot be maintained, the Vesta Charger needs to be monitored to ensure normal function.

Warning: The Vesta Charger must not be used onboard an aircraft.

Warning: Request permission from the ship's crew prior to using your Vesta Charger onboard a ship.

To charge the battery of the OPTIMIZER Smart Mini IPG, perform the following steps:

1. Place the patient in a stationary, comfortable sitting position, ideally reclining at a 45° angle (e.g., on a sofa or armchair).
2. Determine the location of the OPTIMIZER Smart Mini IPG (typically right upper chest area) and then place the flat side of the Vesta charging wand (the side with the four blue rubber screw covers) directly over the OPTIMIZER Smart Mini implant site (over the patient's clothes). To prevent the charging wand from becoming displaced while charging, the charging wand cable may be draped around the patient's neck or the clip on the charging wand cable may be attached to the patient's clothing.
3. Start the charging process by pressing the **Power Button**, holding the button down for 1-2 seconds, and then releasing it.
4. The charging process begins with a display of the IPG Data Download and IPG Data Download Success screens. **See Figures 34 and 35.**




Figure 34: IPG Data Download Screen




Figure 35: IPG Data Download Success Screen

5. After the data download has been completed, the Charging IPG Status screen is displayed by the Vesta Charger. **See Figure 36.**


The Coupling Level icon () at the center of the Charging IPG Status screen will show anywhere from zero to four illuminated bars. Reposition the charging wand until at least 2 bars of the Coupling Level icon are illuminated.



Figure 36: Charging IPG Status Screen

Note: Zero illuminated bars on the Coupling Level icon accompanied by an audible beeping tone indicates poor placement of the charging wand. If the charging wand is not repositioned onto the implant site within 20 seconds, the Vesta Charger will emit 3 long beeping tones, display the Charging IPG Coupling Error screen, and then shut off. If this occurs, press the **Power Button** again to initiate a new charging session.

6. The number of bars on the Charging IPG Battery icon (see icon image on the right) depicts the current charge level of the OPTIMIZER Smart Mini IPG.
7. The Charging IPG Status screen (see **Figure 36**) will continue to be displayed as the OPTIMIZER Smart Mini IPG is being charged.



Note: It is recommended that the patient remain stationary during the charging process. If the charging wand becomes displaced during charging, the Coupling Level icon will show zero illuminated bars and the Vesta Charger will begin to emit an audible beeping tone. If this occurs, please reposition the charging wand until at least 2 bars are illuminated on the Coupling Level icon.

Note: Instruct the patient to try to fully charge their OPTIMIZER Smart Mini IPG during the charging session. Also, inform the patient that charging their implanted device may take longer than one hour if its battery is significantly depleted. If the recharging of the OPTIMIZER Smart Mini IPG cannot be completed in one session, instruct the patient to repeat charging sessions (at least daily) until their implanted device is fully charged.

8. When the battery of the OPTIMIZER Smart Mini IPG is fully charged, the Vesta Charger will emit three short beeping tones and display the IPG Charging Successfully Completed screen (see **Figure 37**). The Vesta Charger will then shut off automatically.




Figure 37: IPG Charging Successfully Completed Screen

9. Detach the charging wand cable clip from the patient's clothing (if necessary), then remove the Vesta charging wand from the patient's implant site and undrape the wand cable from around the patient's neck.
10. Reconnect the AC Adapter to the Vesta Charger as described in Section 1.11.

1.9.1 Early Termination of Charging Session

To terminate a charging session before it has been completed, instruct the patient to press and hold the **Power Button** down for one second and then release it. The Vesta Charger will emit 3 short beeping tones and display the Charge Session Cancellation screen. **See Figure 38.**




Figure 38: Charge Session Cancellation Screen

Alternatively, the patient can remove the charging wand of the Vesta Charger from the implant site, which will cause the Vesta Charger to time out and shut off automatically.

Note: During the charging process, the Vesta Charger monitors the temperature of the OPTIMIZER Smart Mini IPG. To resume charging the OPTIMIZER Smart Mini IPG after terminating a charging session, please wait for approximately 10 minutes before initiating a new charging session to allow the temperature of the implanted device to return to its baseline temperature.

1.10 Charging the OPTIMIZER Smart Mini IPG in Special Charge Mode

If an OPTIMIZER Smart Mini IPG is unable to be charged conventionally because of an alert condition (i.e., Safe Mode), the OPTIMIZER Smart Mini IPG may be charged using the Special Charge Mode.

Note: The OPTIMIZER Smart Mini IPG must be paired with the Vesta Charger before using it to charge the IPG in Special Charge Mode. If necessary, use the instructions in section 1.7 to pair the Vesta Charger with the OPTIMIZER Smart Mini IPG before proceeding.

To charge an OPTIMIZER Smart Mini IPG in Special Charge Mode, perform the following steps:

1. Place the patient in a stationary, comfortable sitting position.
2. Determine the location of the OPTIMIZER Smart Mini IPG (typically right upper chest area) and then place the flat side of the Vesta charging wand (the side with the four blue rubber screw covers) directly over the OPTIMIZER Smart Mini implant site (over the patient's clothes).
3. Place a pairing magnet (or a standard pacemaker magnet) to the left of the **Power Button** on the Vesta Charger. **See Figure 39.**



Figure 39: Pairing Magnet on Vesta Charger

4. Start the charging process by pressing and holding down the **Power Button** (> 5 seconds) until the Vesta Charger emits a single beeping tone, and then releasing it.
5. When the Vesta Charger is used in Special Charge Mode, the charger skips the IPG data download and begins charging the IPG. The Charging IPG Status screen shown in **Figure 40** will be displayed during the charging session.



Figure 40: Charging IPG Status Screen When Charging IPG in Special Charge Mode

1.11 Vesta Charger Placement When Not Being Used for Device Charging

Whenever the Vesta Charger is not being used to charge the OPTIMIZER Smart Mini IPG, instruct the patient to place it in an area frequented by the patient (e.g., bedside table in the bedroom), connected to its AC Adapter, and the AC Adapter plugged into the wall outlet. This will keep the battery of the Vesta Charger fully charged as well as ensure regular communications between the OPTIMIZER Smart Mini IPG and the Vesta Charger.

1.12 Frequency of Charging Sessions

The optimal performance of the rechargeable battery in the OPTIMIZER Smart Mini IPG is only ensured if the battery is fully recharged every week. The day or time chosen to charge the OPTIMIZER Smart Mini IPG is not important, however, it is recommended that the patient not let more than one week pass between the charge sessions.

If the Vesta Charger is not used to perform a charging session on the OPTIMIZER Smart Mini IPG within the time period set by the Optimizer Smart Mini Programmer application, the patient may see the Long Time Without Charging IPG alert screen (see **Figure 41**) displayed by the Vesta Charger.



Figure 41: Long Time Without Charging IPG Alert Screen

If a patient reports seeing this screen displayed by the Vesta Charger, instruct the patient to use their Vesta Charger to charge their OPTIMIZER Smart Mini IPG. If the patient reports that their attempt to charge their OPTIMIZER Smart Mini IPG with their Vesta Charger was unsuccessful, please contact your Impulse Dynamics representative.

If the battery voltage of the OPTIMIZER Smart Mini IPG battery drops below 3.5 V, CCM therapy delivery is automatically suspended. If this occurs, the OPTIMIZER Smart Mini IPG will need to be recharged before it resumes delivering CCM therapy. Once the OPTIMIZER Smart Mini IPG has been recharged, it will automatically resume CCM therapy delivery with its previously programmed parameters.

1.13 Communications

1.13.1 Communications with the OPTIMIZER Smart Mini IPG

The Vesta Charger is configured to communicate with the OPTIMIZER Smart Mini IPG at least once a day. This communication occurs whenever the IPG is within 1.5 m (5 ft) of the Vesta Charger for a few minutes.

If the Vesta Charger and the OPTIMIZER Smart Mini IPG do not communicate within the time period set by the OPTIMIZER Smart Mini Programmer application, the patient may see the Long Time Without Downloading Data From IPG alert screen (see **Figure 42**) displayed by the Vesta Charger.




Figure 42: Long Time Without Downloading Data From IPG Alert Screen

If a patient reports seeing this screen displayed by the Vesta Charger, instruct the patient to attempt to charge their OPTIMIZER Smart Mini IPG with their Vesta Charger. If the patient is able to charge their implanted device successfully, then the alert screen should no longer be displayed by the Vesta Charger. If a patient reports that their attempt to charge their OPTIMIZER Smart Mini IPG with their Vesta Charger was unsuccessful, please contact your Impulse Dynamics representative.

1.14 Call Doctor Alert Codes

In addition to charging the OPTIMIZER Smart Mini IPG, the Vesta Charger is also able to notify the patient of an alert condition in the OPTIMIZER Smart Mini IPG that requires action.

If a detected alert condition is associated with a Direct Action Alert, an alert screen such as Long Time Without Downloading Data from IPG (see **Figure 42**) will be displayed by the Vesta Charger.

If the detected condition is associated with a Call Doctor Alert, the Vesta Charger will display a Call Doctor Alert Code (preceded by a letter denoting the IPG model code) on its screen. The display of a Call Doctor Alert Code (with the exception of code 32) is dependent on whether the specific Patient Alert associated with the Call Doctor Alert Code has been enabled using the Optimizer Smart Mini Programmer application.

Table 5: Call Doctor Alert Codes for the OPTIMIZER Smart Mini IPG

Alert Code	Alert Description	Prevents Charge	Persistent	Auto Refresh
9	IPG Deactivated (see Section 1.14.2.1)	Yes	No	Yes
19	Lead Impedance Change (see Section 1.14.2.2)	No	Yes	Yes
21	CCM Therapy Suspended (see Section 1.14.2.3)	No	No	No
23	Low IPG Battery Voltage (see Section 1.14.2.4)	No	No	Yes
25	CCM Not Sensing/Noise (see Section 1.14.2.5)	No	Yes	Yes
27	Low CCM Therapy Rate (see Section 1.14.2.6)	No	Yes	Yes
31	Charger Failure (see Section 1.14.2.7)	Yes	N/A	N/A
32	IPG has not been paired with Charger (see Section 1.14.2.8)	Yes	N/A	N/A

1.14.1 Call Doctor Alert Code Attributes

Each Call Doctor Alert has the following attributes:

- **Prevents Charge:** An alert that forces the Vesta Charger to terminate the charging process.
- **Persistent:** An alert that will be displayed even if the alert condition that triggered the event is no longer present.
- **Auto Refresh:** An alert that will be displayed again after 24 hours if the alert condition is still present.

1.14.2 Call Doctor Alert Code Definitions

The OPTIMIZER Smart Mini IPG supports the following Call Doctor Alert Codes.

1.14.2.1 Alert Code 9

When Alert Code 9 is displayed, it means that the OPTIMIZER Smart Mini IPG has been deactivated and placed in Safe Mode. If the Vesta Charger displays this Alert Code, please contact your Impulse Dynamics representative.

1.14.2.2 Alert Code 19

When Alert Code 19 is displayed, it means that the OPTIMIZER Smart Mini IPG has detected a significant change in the impedance in one or both ventricular leads. If the Vesta Charger displays this Alert Code, please contact your Impulse Dynamics representative.

1.14.2.3 Alert Code 21

When Alert Code 21 is displayed, it means that CCM therapy in the OPTIMIZER Smart Mini IPG has been suspended. If the Vesta Charger displays this Alert Code, please contact your Impulse Dynamics representative.

1.14.2.4 Alert Code 23

When Alert Code 23 is displayed, it means that the battery voltage level in the OPTIMIZER Smart Mini IPG is less than 3.6 V. If the Vesta Charger displays this Alert Code, please charge the OPTIMIZER Smart Mini IPG battery as soon as possible to prevent CCM therapy from becoming suspended.

1.14.2.5 Alert Code 25

When Alert Code 25 is displayed, it means that the OPTIMIZER Smart Mini IPG has detected that an implanted lead is not sensing or sensing an excessive amount of noise. If the Vesta Charger displays this Alert Code, please contact your Impulse Dynamics representative.

1.14.2.6 Alert Code 27

When Alert Code 27 is displayed, it means that the OPTIMIZER Smart Mini IPG has detected that the amount of CCM therapy delivered is below the alert level programmed into the implanted device by the Optimizer Smart Mini Programmer application. If the Vesta Charger displays this Alert Code, please contact your Impulse Dynamics representative.

1.14.2.7 Alert Code 31

When Alert Code 31 is displayed, it means that the Vesta Charger has detected repeated internal errors during its operation. If the Vesta Charger displays this Alert Code, please contact your Impulse Dynamics representative.

1.14.2.8 Alert Code 32

When Alert Code 32 is displayed, it means that the Vesta Charger has determined that it is attempting to be used on an unrecognized device. If this Alert Code is displayed by the Vesta Charger, please pair the Vesta Charger with the implanted OPTIMIZER Smart Mini IPG and then restart the charging process. If the Vesta Charger still displays this code after it has been successfully paired with the implanted OPTIMIZER Smart Mini IPG, please contact your Impulse Dynamics representative.

1.15 FCE Charger

The FCE Charger allows the field clinical engineer (FCE) or clinical staff to charge a patient's implanted OPTIMIZER Smart Mini IPG in a clinical setting without disrupting its pairing with the patient's assigned Vesta Charger.

Note: FCE Chargers are for clinical use only and are not to be assigned to patients. To obtain an FCE Charger, please contact your Impulse Dynamics representative.

The FCE Charger is intended to be used in a clinical setting to charge a patient's implanted OPTIMIZER Smart Mini IPG whose battery charge level is found to be too low to allow for device interrogation / programming or when the implanted IPG is found to be in DOWN Mode and the battery charge level is too low to allow for an IPG reset.

Using the Unpaired Charge Mode, the FCE Charger can be used in a clinical setting to recharge a deeply discharged OPTIMIZER Smart Mini IPG that cannot be paired and charged with a standard Vesta Charger.

1.15.1 Screens Displayed When FCE Charger is Connected to AC Adapter

1.15.1.1 FCE Charger Self-Charge Status Screen

This screen is displayed whenever the AC Adapter is connected to the FCE Charger. The number of bars shown on the battery icon will vary depending on the current level of charge in the FCE Charger battery (see **Table 2** in section 1.6.1.1).



Figure 43: FCE Charger Self-Charge Status Screen

1.15.1.2 FCE Charger Self-Charge Success Screen

This screen is displayed either when the AC Adapter has successfully completed charging the internal battery of the FCE Charger, when the AC Adapter is connected to the FCE Charger and the battery charge level of the FCE Charger is above 75%, or when the AC Adapter is charging the FCE Charger and the AC Adapter current is less than 50 mA.



Figure 44: FCE Charger Self-Charge Success Screen

1.15.2 Charging the OPTIMIZER Smart Mini IPG using the FCE Charger

To charge an OPTIMIZER Smart Mini IPG using the FCE Charger, perform the following steps:

1. Place the patient in a stationary, comfortable sitting position.
2. Determine the location of the OPTIMIZER Smart Mini IPG (typically right upper chest area) and then place the flat side of the Vesta charging wand (the side with the four blue rubber screw covers) directly over the OPTIMIZER Smart Mini implant site (over the patient's clothes).
3. Start the charging process by pressing the **Power Button**, holding the button down for 1-2 seconds, and then releasing it.

- Using the FCE Charger, the charging process begins by displaying the FCE Charger/IPG Pairing screen as the FCE Charger attempts to pair with the OPTIMIZER Smart Mini IPG. **See Figure 45.**

Note: If the FCE Charger is unable to successfully pair with the patient's implanted OPTIMIZER Smart Mini IPG, reposition the charging wand, and repeat Step 3. If pairing is still unsuccessful, proceed to section 1.15.3.




Figure 45: FCE Charger/IPG Pairing Screen

- When the pairing process has been successfully completed, the FCE Charger will emit 3 short beeping tones and display the Charger/IPG Pairing Success screen. **See Figure 46.**



Figure 46: Charger/IPG Pairing Success Screen

- After the pairing has been successfully completed, the FCE Charger will display the FCE Charging IPG Status screen. **See Figure 47.**


The Coupling Level icon (), at the center of the Charging IPG Status screen will show anywhere from zero to four illuminated bars. Reposition the charging wand until at least 2 bars of the Coupling Level icon are illuminated.



Figure 47: FCE Charging IPG Status Screen

Note: Zero illuminated bars on the Coupling Level icon accompanied by an audible beeping tone indicates poor placement of the charging wand. If the charging wand is not repositioned onto the implant site within 20 seconds, the Vesta Charger will emit 3 long beeping tones, display the Charging IPG Coupling Error screen, and then shut off. If this occurs, press the **Power Button** again to initiate a new charging session.

- The FCE Charging IPG Status screen (see **Figure 47**) will continue to be displayed as the OPTIMIZER Smart Mini IPG is being charged.

Note: It is recommended that the patient remain stationary during the charging process.

Note: If the intended use of the FCE Charger is only to charge the battery of the OPTIMIZER Smart Mini IPG enough to allow for the interrogation / programming of the device, the charging session may be terminated once the IPG Battery icon (see icon on the right in **Figure 47**) displays 2 bars, with the last one flashing. To terminate the charge session, press the **Power Button**, hold it down for 1-2 seconds, and then release it.

- When the battery of the OPTIMIZER Smart Mini IPG is fully charged, the Vesta Charger will emit three short beeping tones and display the IPG Charging Successfully Completed screen (see **Figure 48**). The Vesta Charger will then shut off automatically.




Figure 48: IPG Charging Successfully Completed Screen

- Remove the Vesta charging wand from the patient's implant site and undrape the wand cable from around the patient's neck.

1.15.2.1 Charging the OPTIMIZER Smart Mini IPG in Safe Mode

When charging an OPTIMIZER Smart Mini IPG that is in Safe Mode, the FCE Charging IPG Status screen shown in **Figure 49** will be displayed during the charging session.



Figure 49: FCE Charging IPG Status Screen When Charging IPG in Safe Mode

1.15.3 Charging the OPTIMIZER Smart Mini IPG in Unpaired Charge Mode

To charge an OPTIMIZER Smart Mini IPG using the FCE Charger in Unpaired Charge Mode, perform the following steps:

- Place the patient in a stationary, comfortable sitting position.
- Determine the location of the OPTIMIZER Smart Mini IPG (typically right upper chest area) and then place the flat side of the Vesta charging wand (the side with the four blue rubber screw covers) directly over the OPTIMIZER Smart Mini implant site (over the patient's clothes).
- Place a pairing magnet (or a standard pacemaker magnet) to the left of the **Power Button** on the Vesta Charger. **See Figure 50.**



Figure 50: Pairing Magnet on FCE Charger

4. Start the unpaired charging process by pressing and holding down the **Power Button** (> 10 seconds) until the Vesta Charger automatically starts the charging process, and then releasing it.
5. When using the FCE Charger in Unpaired Charge Mode, the charging process begins by displaying the FCE Charging IPG Status screen. **See Figure 51.**



Figure 51: FCE Charging IPG Status Screen in Unpaired Charge Mode

6. If the FCE Charger is able to reestablish communications with the OPTIMIZER Smart Mini IPG within 5 minutes of charging in Unpaired Charge Mode, the FCE Charger will emit three short beeping tones, display the Charge Session Cancellation screen (see **Figure 52**), and then automatically terminate the Unpaired Charge Mode charging session. When this occurs, remove the Pairing Magnet from atop the FCE Charger and then proceed Step 3 in section 1.15.2.

Note: If the FCE Charger is not able to reestablish communications with the OPTIMIZER Smart Mini IPG after 5 minutes of charging in Unpaired Charge Mode, it will emit three short beeping tones, display the Charge Session Cancellation screen (see **Figure 52**), and then terminate the Unpaired Charge Mode charging session. If this occurs, repeat Step 4 to initiate another charge session in Unpaired Charge Mode.




Figure 52: Charge Session Cancellation Screen

1.16 Cleaning

Warning: Always unplug the AC Adapter from the Vesta Charger prior to cleaning.

The exterior surface of the Vesta Charger should only be cleaned with disinfectant wipes as needed.

Caution: DO NOT use solvents or cleaning cloths impregnated with chemical cleaning agents.

Warning: DO NOT attempt to clean the electrical connector of your Vesta Charger.

Warning: DO NOT submerge any part of the Vesta Charger in water. Damage to the unit may result. The Vesta Charger has limited protection against the ingress of water or humidity (ingress protection rating IP22).

Warning: DO NOT sterilize any part of the Vesta Charger because any such attempt could severely damage the equipment.

1.17 Maintenance

The Vesta Charger does not contain any user-serviceable parts. If the Vesta Charger is not operational, please contact your Impulse Dynamics representative to obtain a replacement charger.

Warning: No modification of this equipment is allowed.

The battery inside the Vesta Charger is expected to have a service life of at least 5 years. If the Vesta Charger cannot fully charge an OPTIMIZER Smart Mini IPG after the charger battery has been fully charged, please contact your Impulse Dynamics representative to obtain a replacement charger.

1.18 Storage and Handling

The Vesta Charger System is designed to maintain functionality after it has been exposed to the following environmental extremes:

- **Ambient Temperature:** -20°C to +60°C (-4°F to 140°F)
- **Relative Humidity:** 10% to 100% (with or without condensation)
- **Atmospheric Pressure:** 50 kPa to 156 kPa (14.81 inHg to 46.20 inHg)

The Vesta Charger System should not be exposed to excessively hot or cold storage conditions. Patients should be instructed not to leave the Vesta Charger System in their car or outdoors for extended periods of time. Temperature extremes, particularly high heat can damage the sensitive electronics of the Vesta Charger System.

For proper operation, the Vesta Charger should be used only under the following environmental conditions:

- **Ambient Temperature:** 10°C to 27°C (50°F and 81°F)
- **Relative Humidity:** 20% to 75%
- **Atmospheric Pressure:** 70 kPa to 106 kPa (20.73 inHg to 31.39 inHg)

Note: When not being used to charge the OPTIMIZER Smart Mini IPG, the Vesta Charger should always be connected to its AC Adapter, and the AC Adapter plugged into the wall outlet

1.19 Disposal

If the Vesta Charger is no longer needed by the patient and is returned, please notify your Impulse Dynamics representative of its return.

Warning: **DO NOT** discard the Vesta Charger in the trash bin. The Vesta Charger contains Lithium batteries as well as non-RoHS components. If disposal of the Vesta Charger is necessary, properly dispose of the Vesta Charger in accordance with local regulations governing the disposal of such material.


THIS PAGE INTENTIONALLY LEFT BLANK

APPENDIX I

Electromagnetic Immunity

Electromagnetic Immunity of the Vesta Charger

GUIDELINES AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY OF THE VESTA CHARGER			
<p>Essential Performance of the Vesta Charger:</p> <ul style="list-style-type: none"> • The Vesta Charger shall not charge any IPG inappropriately. • The Vesta Charger shall only charge a paired IPG appropriately. • The patient shall be made aware of inappropriate charging either by an explicit message, or by the absence of an expected message from the Vesta Charger. • If essential performance is lost due to electromagnetic disturbances, the Vesta Charger shall not be able to charge any IPG. 			
<p>The Vesta Charger, part of the OPTIMIZER Smart Mini System, is intended for use in an electromagnetic environment as specified below. The customer or user of the Vesta Charger must ensure that it is used within the specified environment.</p> <p>The test levels follow FDA recommendations for the home environment per “Design Considerations for Devices Intended for Home Use - Guidance for Industry and Food and Drug Administration Staff”, November 24, 2014</p>			
Immunity test	IEC 60601-1-2:2014 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic discharge as defined in IEC 61000-4-2	Contact Discharge: ± 8 kV Air Discharge: ± 2 kV, ± 4 kV, ± 8 kV, and ± 15 kV	Contact Discharge: ± 8 kV Air Discharge: ± 2 kV, ± 4 kV, ± 8 kV, and ± 15 kV	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, relative humidity should be 30% or greater.
Electrical fast transient / burst as defined in IEC 61000-4-4	± 2 kV for mains power supply ± 1 kV for in-/output lines	± 2 kV for mains power supply ± 1 kV for in-/output lines	Mains power quality should be that of a typical home healthcare, business, or hospital environment. Do not operate motors or other noisy electrical equipment on the same mains circuit as the Vesta Charger.
AC line voltage surges as defined in IEC 61000-4-5	± 2 kV Common Mode ± 1 kV Differential Mode 1.2/50 µs	± 2 kV Common Mode ± 1 kV Differential Mode 1.2/50 µs	Mains power quality should be that of a typical home healthcare, business, or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines as defined in IEC 61000-4-11	0%, 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0%, 1 cycle 70%, 25 cycles 0%, 250 cycles	0%, 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0%, 1 cycle 70%, 25 cycles 0%, 250 cycles	Mains power quality should be that of a typical home healthcare, business, or hospital environment. Note: If the user of the Vesta Charger requires uninterrupted operation during power mains interruptions, it is recommended to power the Vesta Charger from an uninterruptible power supply.

Power line frequency magnetic fields (50/60 Hz) as defined in IEC 61000-4-8	30 A/m	30 A/m	Power line frequency magnetic fields (50/60 Hz) should be at levels expected in a typical home healthcare, business, or hospital environment.
Conducted RF as defined in IEC 61000-4-6:2013	3 V r.m.s outside industrial, scientific, and medical (ISM) and amateur radio bands between 0.15 MHz and 80 MHz, 6 V r.m.s. in ISM and amateur radio bands between 0.15 MHz and 80 MHz	3 V r.m.s outside industrial, scientific, and medical (ISM) and amateur radio bands between 0.15 MHz and 80 MHz, 6 V r.m.s. in ISM and amateur radio bands between 0.15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Radiated RF as defined in IEC 61000-4-3: 2006 +A1: 2007 +A2: 2010	10 V/m: 80 MHz to 2.7 GHz, 80% 1kHz AM	10 V/m: 80 MHz to 2.7 GHz, 80% 1kHz AM	$d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P}$ 80 MHz to 800 MHz $d = 2.33\sqrt{P}$ 800 MHz to 2.5 GHz Where "P" is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and "d" is the recommended separation distance in meters (m).
Proximity fields from RF communications equipment as defined in IEC 61000-4-3	Various per table 9	Various per table 9	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, "a" should be less than the compliance level in each frequency range "b". Interference may occur in the vicinity of equipment marked with the following symbol: 
Notes: a - Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be theoretically predicted with accuracy. An electromagnetic site survey should be taken into consideration to assess the electromagnetic environment due to fixed RF transmitters. If the measured field strength in the location where the Vesta Charger is used exceeds the applicable RF compliance level above, the Vesta Charger should be monitored to ensure normal operation. If an abnormal function is observed, additional measures may be necessary, such as relocating the Vesta Charger. b - For frequencies in the range of 150 kHz to 80 MHz, the field strength should be less than 3 V/m.			

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the Vesta Charger

Recommended separation distances between portable and mobile RF communications equipment and the Vesta Charger			
The Vesta Charger should be used in an electromagnetic environment with limited radiated RF noise. The customer or user of the Vesta Charger can help prevent electromagnetic interference by maintaining the minimum distance between portable and mobile RF communications equipment (transmitters) and the Vesta Charger recommended below, which is determined by the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance broken down by transmitter frequency(m)		
	150 kHz to 80 MHz¹ d = 1.17√P	80 MHz to 800 MHz¹ d = 1.17√P	800 MHz to 2.5 GHz d = 2.33√P
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.75
1	1.17	1.17	2.33
10	3.70	3.70	7.36
100	11.70	11.70	23.30
For transmitters with a maximum rated output power not listed above, the recommended separation distance “d” in meters (m) can be estimated by using the equation applicable to the frequency of the transmitter, where “P” is the maximum rated output power of the transmitter in watts (W) specified by the transmitter manufacturer.			
¹ At 80 MHz and 800 MHz, the higher frequency range applies.			
Note: These guidelines may not apply to all settings. Electromagnetic propagation is affected by absorption and reflection from buildings, objects, and people.			

Electromagnetic Emissions

Electromagnetic Emissions from the Vesta Charger

The Vesta Charger must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.

Warning: The Vesta Charger must not be used onboard an aircraft.

Warning: Permission must be requested from a ship's crew prior to using the Vesta Charger onboard a ship.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Vesta Charger, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

ETSI EN 300 330

GUIDELINES AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS OF THE VESTA CHARGER PURSUANT TO:		
ETSI EN 300 330 V2.1.1 – Short Range Devices (SRD); Radio equipment in the frequency range 9kHz to 25MHz and inductive loop systems in the frequency range 9kHz to 30MHz; Harmonised Standard covering the essential requirements of article 3.2 of Directive 2014/53/EU		
The Vesta Charger, part of the OPTIMIZER Smart Mini System, is intended for use in an electromagnetic environment as specified below. The customer or user of the Vesta Charger must ensure that it is used within the specified environment.		
Emissions Test	Compliance	Electromagnetic environment – guidelines
Permitted range of operating frequencies	Complies with clause 4.3.2.3	The Vesta Charger must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
Modulation bandwidth	Complies with clause 4.3.3.3	
Radiated H Field	Complies with clause 4.3.4.3	
Transmitter Spurious Emissions below 30MHz – Operating and Stand-By Mode	Complies with clause 4.3.8.3	
Transmitter Spurious Emissions 30-1000MHz – Operating and Stand-By Mode	Complies with clause 4.3.9.3	
Receiver Spurious Emissions up to 1000MHz	Complies with clause 4.4.2.3	

ETSI EN 301 839

<p>GUIDELINES AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS OF THE VESTA CHARGER PURSUANT TO:</p> <p>ETSI EN 301 839 V2.1.1 – Ultra Low Power Active Medical Implants (ULP-AMI) and associated Peripherals (ULP-AMI-P) operating in the frequency range 402 MHz to 405 MHz; Harmonised Standard covering the essential requirements of article 3.2 of the Directive 2014/53/EU</p>		
<p>The Vesta Charger, part of the OPTIMIZER Smart Mini System, is intended for use in an electromagnetic environment as specified below. The customer or user of the Vesta Charger must ensure that it is used within the specified environment.</p>		
Emissions Test	Compliance	Electromagnetic environment – guidelines
Frequency Error	Complies with clause 4.2.1.1	The Vesta Charger must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
Emission Bandwidth	Complies with clause 4.2.1.2	
EIRP	Complies with clause 4.2.1.3	
Transmitter Spurious Emissions (30 MHz to 6 GHz)	Complies with clause 5.3.4	
Frequency Stability Under Low Voltage Conditions	Complies with clause 4.2.1.5	
Receiver Spurious Emissions	Complies with clause 4.2.2.1	
Spectrum Access	Complies with clause 4.2.3.1	
Receiver Blocking	Complies with clause 4.2.3.2	

ETSI EN 301 489-1 and ETSI EN 301 489-27

<p>GUIDELINES AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS OF THE VESTA CHARGER PURSUANT TO:</p> <p>ETSI EN 301 489-1 V2.2.3 – ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements; Harmonised Standard for ElectroMagnetic Compatibility</p> <p>ETSI EN 301 489-27 – ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 27: Specific conditions for Ultra Low Power Active Medical Implants (ULP-AMI) and related peripheral devices (ULP-AMI-P) operating in the 402 MHz to 405 MHz bands; Harmonised Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU</p>			
<p>The Vesta Charger, part of the OPTIMIZER Smart Mini System, is intended for use in an electromagnetic environment as specified below. The customer or user of the Vesta Charger must ensure that it is used within the specified environment.</p>			
<p>There may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.</p>			
Emissions Test	Basic Standard	Compliance	Electromagnetic environment – guidelines
Radiated Emissions	EN 55032	N/A – covered by relevant radio standards	The Vesta Charger must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
Conducted Emissions	EN 55032	Pass	
Harmonic Current Emissions	IEC 61000-3-2	Pass	
Voltage Fluctuations	IEC 6100-3-3	Pass	

IEC 60601-1-2 2014

GUIDELINES AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS OF THE VESTA CHARGER PURSUANT TO:

IEC 60601-1-2 2014, Edition 4.0 – Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances – Requirements and tests

The Vesta Charger, part of the OPTIMIZER Smart Mini System, is intended for use in an electromagnetic environment as specified below. The customer or user of the Vesta Charger must ensure that it is used within the specified environment.

There may be potential difficulties in ensuring electromagnetic compatibility in other environments, due to conducted as well as radiated disturbances.

Emissions Test	Standard/Section	Compliance	Electromagnetic environment – guidelines
Radiated Emissions 30-1000MHz	CISPR11, Section 6, Table 11 (Class B, Group 2)	Group 2, Class B	The Vesta Charger must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
Conducted Emissions 0.15-30MHz, 230V 50Hz and 120V, 60Hz	CISPR11, Section 6, Table 6 (Class B, Group 2)	Group 2, Class B	
AC Harmonic Emissions	IEC 61000-3-2	Class A	
Voltage Fluctuations	IEC 61000-3-3	Pass	

APPENDIX II

Wireless Technology

RF wireless technology is used to transcutaneously transmit energy from the Vesta Charger to recharge the OPTIMIZER Smart Mini IPG at the 13.56 MHz ISM frequency. The transmission range is specified at a maximum of 4 cm (1.5 in) between the Charger's coil and the IPG's receiving coil. Control over the recharge process, as well as the communications of alert messages from the IPG to the Charger take place over an encrypted MICS channel.

Vesta Charger Wireless Nominal Specifications


Characteristic	Nominal
MICS MedRadio	
Frequency Band	402 – 405 MHz Medical Implant Communication Service (MICS) Medical Device Radio Communication Service (MedRadio)
Bandwidth	240 kHz
Modulation	FSK
Radiated Power	-20.6 dBm EIRP
Range	0 to at least 1.5 m
Transcutaneous Energy Transfer	
Frequency Band	13.56 MHz Industrial, Scientific, and Medical radio band (ISM)
Modulation	Amplitude (slow to optimize coupling)
Radiated Power	< 0.6 W reactive near-field
Range	5 mm to 40 mm
Recharge Channel Communication	
Frequency Band	13.56 MHz \pm 9.2 ppm Industrial, Scientific, and Medical radio band (ISM)
Bandwidth	<0.014 MHz
Modulation	PPM
Radiated Power	-6.93 dBm EIRP
Range	5 mm to 40 mm

Quality of Service (QoS) for Wireless Technology


QoS for Communications between the Vesta Charger and the OPTIMIZER Smart Mini IPG


MedRadio in the MICS sub-band (402 to 405 MHz) wireless technology enables communication between the OPTIMIZER Smart Mini IPG and the Vesta Charger. The requirements for the Quality of Service (QoS) vary depending on the use environment (operating room, recovery room, clinic, and home environment).

The Vesta Charger will begin by displaying the IPG Data Download and IPG Data Download Success screens:



After the data download has been completed, the Charging IPG Status screen is displayed by the Vesta Charger:




The Coupling Level icon (), whose number of illuminated bars is proportional to the proximity of the charging wand to the implanted OPTIMIZER Smart Mini IPG, is indicative of the Quality of Service (QoS) for the transcutaneous energy transmission wireless link. The charging wand should be repositioned until at least 2 bars of the Coupling Level icon are illuminated, indicating sufficient QoS for charging the OPTIMIZER Smart Mini IPG.

One illuminated bar indicates degraded QoS which may require a longer charging time. Zero illuminated bars on the Coupling Level icon accompanied by an audible beeping tone indicates poor placement of the charging wand. If the charging wand is not repositioned onto the implant site within 20 seconds, the Vesta Charger will emit 3 long beeping tones, display the Charging IPG Coupling Error screen, and then shut off.

Besides charging the OPTIMIZER Smart Mini, the Vesta Charger also serves as a way of messaging the patient about alerts and other conditions. The Vesta Charger is configured to communicate with the OPTIMIZER Smart Mini IPG at least once a day. This communication occurs whenever the IPG is within 1.5 m (5 ft) of the Vesta Charger for a few minutes.

If the Vesta Charger and the OPTIMIZER Smart Mini IPG do not communicate within a programmable time period, the patient may see the “Long Time Without Downloading Data From IPG” alert screen displayed by the Vesta Charger:



In this case, instruct the patient to attempt to charge their OPTIMIZER Smart Mini IPG with their Vesta Charger. If the patient is able to charge their implanted device successfully, then the alert screen should no longer be displayed by the Vesta Charger. If the attempt to charge the OPTIMIZER Smart Mini IPG with the Vesta Charger is unsuccessful, the Impulse Dynamics representative should be contacted.

Wireless Security Measures

Wireless Security Measures in Communications between OPTIMIZER Smart Mini IPG and Vesta Charger

Pairing of the Vesta Charger with the OPTIMIZER Smart Mini IPG ensures that the communication and charging information received by the Vesta Charger is securely encrypted and unique to a specific implanted device.

During the pairing process, the Vesta Charger uses short-range communication to search for a device to pair with and creates an encryption key once a compatible device model has been found. This encryption key is stored and used by the Vesta Charger for all its subsequent communications sessions with the paired device.

Wireless signals are secured through device system design that includes the following:

- Pairing of a Vesta Charger and an OPTIMIZER Smart Mini IPG requires the placing of a pairing magnet on the Vesta Charger and locating the Charge Wand within 4 cm (1.5 in) of the OPTIMIZER Smart Mini IPG. The 13.56 MHz short-range channel is used as part of a proprietary process to pair the devices and exchange encryption keys.
- The OPTIMIZER Smart Mini IPG and the Vesta Charger encrypt their wireless communications using encryption keys that are generated during the pairing process.
- Only one Vesta Charger can be paired with the IPG at any one time.

Troubleshooting for Wireless Coexistence Issues

Troubleshooting Wireless Connection between OPTIMIZER Smart Mini IPG and Vesta Charger

If you experience issues with establishing a wireless connection between the OPTIMIZER Smart Mini IPG and the Vesta Charger, try the following:

- Whenever the Vesta Charger is not being used to charge the OPTIMIZER Smart Mini IPG, place it in an area that is frequented by the patient (e.g., bedside table in the bedroom), connected to its AC Adapter, and the AC Adapter plugged into the wall outlet. This will ensure regular communications between the OPTIMIZER Smart Mini IPG and the Vesta Charger.
- Remain stationary during the charging or data transfer process.
- Decrease the distance between the devices.
- Move the devices so they share line of sight.
- Move the devices away from other devices that may be causing interference.
- Do not operate other wireless devices (i.e., programmers for other devices, laptop, tablet, mobile phone, or cordless phone) at the same time.
- Wait a few minutes and try connecting again.

Note: Wireless communications equipment, such as wireless home network devices, mobile and cordless telephones, and tablets, could affect the quality of the wireless connection.